


Size= 134X27.5 mm
 Printing area=121.6X18.5 mm
 Colour of print = Blue (Pantone 2728cvc)
 Pitch=29mm, Cavity=6
 Pitch=32mm, Cavity=5
 Pitch=28.6mm, Cavity=14

HEALFLON[®]
 I.V. Cannula with Catheter &
 Injection Valve

PEEL OPEN



Harsoria

Manufactured By:
HARSORIA HEALTHCARE PVT. LTD.
 110 - 111, PHASE - IV, UDYOG VIHAR,
 GURGAON-122015, HARYANA, INDIA
 Customer Care No.: +91 124 4523400
M.L.No.: MFG/MD/2019/000135

EU REP Mdi Europa GmbH AW/B_01, Rev. 13
 Langenhagener Str. 71
 30855 Langenhagen
 GERMANY - Email: info@mdi-europa.com

UDI



Non-pyrogenic
 STERILE
 MD
 Non-Toxic/No
 Hazardous
 Substances
 DEHP
 LATEX
 CE 0123

(01)XXXXXXXXXXXXXX
 (17)YYMMDD
 (10)XXXXXX/XXXX

REF LOT

Ø:

FLOW:

FILE NAME/TYPE	HEALFLON (TF_01) / FLON TYPE		 Harsoria HARSORIA HEALTHCARE PVT. LTD. CMYK		
DESIGNER	BALBEER SINGH		PANTONE COLORS USED BLUE 2728		
MINIMUM FONT TYPE AND SIZE	Current Size 3.5	Standard Size 4.5			
	Arial	Time New Roman/Arial			
DATE STARTED	12/09/2025				
DATE COMPLETE	01/10/2025				
FILE LOCATION	Z:\CANNULA\FLON TYPE\HEALFLON (TF-01)\B_01, Rev. 13.cdr		SIZE (LXH) -134X27.5mm		



Size= 114x27.5 mm
 Printing area= 104x18 mm
 Color of print= Blue (Pantone2728 cvc)
 Pitch= 29mm, Cavity=6
 Pitch= 32mm, Cavity=5
 Pitch= 28.6mm, Cavity=14



FILE NAME/TYPE	HEALFLON (TF_01) / FLON TYPE		 HARSORIA HEALTHCARE PVT. LTD. CMYK		
DESIGNER	BALBEER SINGH		PANTONE BLUE 2728 COLORS USED		
MINIMUM FONT TYPE AND SIZE	Current Size	Standard Size			
	3	4.5			
	Arial	Time New Roman/Arial			
DATE STARTED	12/09/2025				
DATE COMPLETE	01/10/2025				
FILE LOCATION	D:\Design-Artwork\CANNULA\FLON TYPE\HEALFLON (TF-01)\OK\HEALFLON SMOLL (TF_01a)\HEALFLON\B_01a, Rev. 12.cdr		SIZE (LXH) - 114X27.5mm		

FILE NAME/TYPE	HEALFLON (TF_01) / FLON TYPE		HARSORIA HEALTHCARE PVT. LTD. CMYK	
DESIGNER	BALBEER SINGH		PANTONE BLACK 2728 CVC	
MINIMUM FONT TYPE AND SIZE	Current Size	Standard Size	COLORS	
	4.5	4	USED	
DATE STARTED	12/09/2025		SIZE (HxL) - 168X136X150 mm	
DATE COMPLETE	01/10/2025			
FILE LOCATION	Z:\CANNULA\FLON TYPE\HEALFLON (TF-01)OK\IC_01, Rev.19			

Colour of print = Blue (Pantone 2728CVC) & Black
 Colour Scheme = As per design
 Size = 168X136X150 mm
 Print Instruction sheet inside top flap



BACK SIDE

Without Varnish This area

PACKED
 CHECKED
 TO BE PRINT ON
 BOTTOM FLAP IN
 BLACK COLOR

PACKED
 CHECKED
 TO BE PRINT ON
 BOTTOM FLAP IN
 BLACK COLOR

HEALFLON®

I.V. Cannula with Catheter & Injection Valve

MATERIAL USED :

- PP, POM, HDPE, LDPE, Silicon Rubber & Stainless Steel
- Catheter: PTFE/FEP/PUR/ETFE (see the product ref. code)

COATING MATERIAL :

- Silicone Dispersion

DESCRIPTION :

Catheter gauge size & length are identified on the product packaging. The color of the catheter's plastic hub also indicates the gauge size of catheter.

GAUGE: 14G 16G 18G 20G 22G 24G

COLOR: Orange Gray Green Pink Blue Yellow

The materials used to manufacture this device do not contain natural rubber latex or PVC derivatives.

INTENDED PURPOSE:

"Intravenous Cannula is a device, whereby fluids can be introduced into the human circulatory system during infusion therapy."

INDICATIONS :

- Infusion of I.V. Solutions (To maintain hydration and/or correct dehydration if patient is unable to take sufficient volume of oral fluid)
- Intermittent intravenous Drugs administration

CONTRAINDICATIONS:

- Product should not be used in patients with known hypersensitivity to any of the material used including coating materials
- Product should not be used for Administration of high viscous fluids
- Product should not be used for Large volume blood transfusion
- There is no known reactions between the catheter & Magnetic Resonance Imaging (MRI)

INSTRUCTIONS FOR USE:

- The device is manufactured & tested in accordance with the international standard "Over needle Peripheral Catheters EN ISO 10555-1 & 5"
- The connecting port of this device is in compliance to ISO 80369-7
- Carefully select and aseptically prepare the site
- Select suitable size of I.V. Cannula & inspect visually to ascertain that package is intact
- Remove cannula from sterile packing and remove its needle cover
- Grip the cannula from injection port cap & projection provided on hub
- Perform venipuncture & check for flash back of blood in flash back chamber
- Advance the catheter into vein and simultaneously withdraw the needle
- NEVER TRY TO REINSERT THE PARTIALLY OR COMPLETELY WITHDRAWN NEEDLE
- Withdraw the needle completely while pressing the vein just after the tip of catheter into the vein & discard the needle in an appropriate container
- Connect to the I.V. infusion set line
- Cover the puncture site with sterile dressing
- Drugs can be injected with the help of syringe without needle through integrated injection port after removing the port cap. Close the injection port cap after use
- Perform routine monitoring & venipuncture site maintenance according to medical norms
- Based on Clinical Evaluation Report following recommendations are made for use of the device:
- Upper limb placements are preferable to lower limbs
- For blood sampling, it is recommended that larger gauge catheters be used than for infusion
- Use of specialized infusion teams for insertion and monitoring of IV catheters has been shown to lead to better patient outcomes
- In case of any adverse events occurred during the usage of device it must be reported to manufacturer and the national competent authority

INTENDED PATIENT OF USE:

The device can be used in all age group patients.

12G ,13G & 14G to 22G for adults, 24G for child, 26G for neonates.

Note: •Consult doctor before use of this device in pregnant & Breast feeding women.

•Consult doctor for appropriate size / type of the device or depending upon the patient's condition.

UNDESIRABLE SIDE EFFECTS MAY BE CAUSED BY USED OF THE DEVICE :

- Allergy • Inflammation • Swelling • Infection

TERM OF USE : •Recommended Maximum Duration of Use: 96 hours

WARNINGS:

- The use of this product is restricted to a qualified doctor or a trained paramedic.
 - The device is for Single Use Only. Discard after use.
 - Read instructions before use. The product should be used according to the instructions for use.
 - Do not put the Device to Use after the Use by Date or Date of Expiry
 - If there is any change in expected performance of the device or in case of any malfunction the device should be immediately removed & sent back to supplier for analysis
 - For Known/reported adverse events associated with use of this device, refer to the Clinical Evaluation Report .
 - HARSORIA HEALTHCARE DISCLAIMS ANY RESPONSIBILITY FOR POSSIBLE CONSEQUENCES RESULTING FROM IMPROPER USE
 - Any device that is connected to this product must comply with ISO 80369-7 in order to achieve the intended performance of this product & to avoid leakage in the connection
 - The product should not be reprocessed as reprocessing may compromise the integrity of the device and /or lead to device failure.
 - Visually inspect and carefully check the product and packaging before use. Improper transport and handling may cause structural and/or functional damage to device or packaging
 - Do not clean or re-sterilize the device as this compromise the device performance (functionality), integrity and biological features. This may increase the risk of cross contamination or severe deterioration in health and safety of patients.
 - The product is guaranteed sterile & non-pyrogenic if the package has not been damaged or unintentionally opened before use. If package is damaged or unintentionally opened before use, discard as per country's disposal regulation.
 - Do not expose to heat or direct sunlight
 - The product should be used immediately after opening the packaging.
 - Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or infections
- Disposal/ Discard:**
- Dispose off/ Discard the used device in accordance with your Country's Healthcare and Safety Regulations.
 - Used Sharps/ needles should be disposed in sharps disposal container.
 - Used product/ catheter should be disposed in sanitary container or as hospital protocol to prevent possible contamination and cross infection.

INSTRUCTIONS FOR USE

REF

Product Ref. No.

LOT

Batch Number



Date of Manufacture



Use-by-date



Do Not Reuse



Sterilized using Ethylene Oxide and Single Sterile Barrier System



Caution



Consult Instructions for Use



Non-Pyrogenic



Do Not Use if package is Damaged and consult instructions for use



Do not resterilize



Keep away from sunlight



DEHP Free



Latex Free



Temperature limit

MD

Medical Device

UDI

Unique Device Identifier



Keep Dry



Non-Toxic/No Hazardous Substances



Importer



Distributer



Manufactured By:
Harsoria Healthcare Pvt. Ltd.
110 -111, Phase-IV, Udyog Vihar,
Gurgaon - 122 015, Haryana, INDIA.
Email: info@harsoria.com
Website: www.harsoria.com
Customer Care No.: +91 124 4523400

CE 0123

EU REP

Mdi Europa GmbH
Langenhagener Str. 71
30855 Langenhagen - GERMANY
Email: info@mdi-europa.com
Rev. Date :12.09.2025
Ref. Code: AW/IC_01,Rev.19

Without Varnish This area


HEALFLON®



TO BE PRINT ON
BOTTOM FLAP IN
BLACK COLOR

TO BE PRINT ON
BOTTOM FLAP IN
BLACK INK

BACK SIDE



CUTT

HEALFLON®

I.V. Cannula with Catheter & Injection Valve

MATERIAL USED :

- PP, POM, HDPE, LDPE, Silicon Rubber & Stainless Steel
- Catheter: PTFE/FEP/PUR/ETFE (see the product ref. code)

COATING MATERIAL :

- Silicone Dispersion

DESCRIPTION :

Catheter gauge size & length are identified on the product packaging. The color of the catheter's plastic hub also indicates the gauge size of catheter.

GAUGE: 14G 16G 18G 20G 22G 24G

COLOR: Orange Gray Green Pink Blue Yellow

The materials used to manufacture this device do not contain natural rubber latex or PVC derivatives.

INTENDED PURPOSE:

"Intravenous Cannula is a device, whereby fluids can be introduced into the human circulatory system during infusion therapy."

INDICATIONS :

- Infusion of I.V. Solutions (To maintain hydration and/or correct dehydration if patient is unable to take sufficient volume of oral fluid)
- Intermittent intravenous Drugs administration

CONTRAINDICATIONS:

- Product should not be used in patients with known hypersensitivity to any of the material used including coating materials
- Product should not be used for Administration of high viscous fluids
- Product should not be used for Large volume blood transfusion
- There is no known reactions between the catheter & Magnetic Resonance Imaging (MRI)

INSTRUCTIONS FOR USE:

- The device is manufactured & tested in accordance with the international standard "Over needle Peripheral Catheters EN ISO 10555-1 & 5"
- The connecting port of this device is in compliance to ISO 80369-7
- Carefully select and aseptically prepare the site
- Select suitable size of I.V. Cannula & inspect visually to ascertain that package is intact
- Remove cannula from sterile packing and remove its needle cover
- Grip the cannula from injection port cap & projection provided on hub
- Perform venipuncture & check for flash back of blood in flash back chamber
- Advance the catheter into vein and simultaneously withdraw the needle
- NEVER TRY TO REINSERT THE PARTIALLY OR COMPLETELY WITHDRAWN NEEDLE
- Withdraw the needle completely while pressing the vein just after the tip of catheter into the vein & discard the needle in an appropriate container
- Connect to the I.V. infusion set line
- Cover the puncture site with sterile dressing
- Drugs can be injected with the help of syringe without needle through integrated injection port after removing the port cap.
- Close the injection port cap after use
- Perform routine monitoring & venipuncture site maintenance according to medical norms
- Based on Clinical Evaluation Report following recommendations are made for use of the device:
- Upper limb placements are preferable to lower limbs
- For blood sampling, it is recommended that larger gauge catheters be used than for infusion
- Use of specialized infusion teams for insertion and monitoring of IV catheters has been shown to lead to better patient outcomes
- In case of any adverse events occurred during the usage of device it must be reported to manufacturer and the national competent authority

INTENDED PATIENT OF USE:

The device can be used in all age group patients.
12G ,13G & 14G to 22G for adults, 24G for child, 26G for neonates.

Note: •Consult doctor before use of this device in pregnant & Breast feeding women.

•Consult doctor for appropriate size / type of the device or depending upon the patient's condition.

UNDESIRABLE SIDE EFFECTS MAY BE CAUSED BY USED OF THE DEVICE :

- Allergy • Inflammation • Swelling • Infection

TERM OF USE : •Recommended Maximum Duration of Use: 96 hours

pital protocol to prevent possible contamination and cross infection.

WARNINGS:

- The use of this product is restricted to a qualified doctor or a trained paramedic.
- The device is for Single Use Only. Discard after use.
- Read instructions before use. The product should be used according to the instructions for use.
- Do not put the Device to Use after the Use by Date or Date of Expiry
- If there is any change in expected performance of the device or in case of any malfunction the device should be immediately removed & sent back to supplier for analysis
- For Known/reported adverse events associated with use of this device, refer to the Clinical Evaluation Report .
- HARSORIA HEALTHCARE DISCLAIMS ANY RESPONSIBILITY FOR POSSIBLE CONSEQUENCES RESULTING FROM IMPROPER USE
- Any device that is connected to this product must comply with ISO 80369-7 in order to achieve the intended performance of this product & to avoid leakage in the connection
- The product should not be reprocessed as reprocessing may compromise the integrity of the device and /or lead to device failure.
- Visually inspect and carefully check the product and packaging before use. Improper transport and handling may cause structural and/or functional damage to device or packaging
- Do not clean or re-sterilize the device as this compromise the device performance (functionality), integrity and biological features.
- This may increase the risk of cross contamination or severe deterioration in health and safety of patients.
- The product is guaranteed sterile & non-pyrogenic if the package has not been damaged or unintentionally opened before use.
- If package is damaged or unintentionally opened before use, discard as per country's disposal regulation.
- Do not expose to heat or direct sunlight
- The product should be used immediately after opening the packaging.
- Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or infections

Disposal/ Discard:

- Dispose off/ Discard the used device in accordance with your Country's Healthcare and Safety Regulations.
- Used Sharps/ needles should be disposed in sharps disposal container.
- Used product/ catheter should be disposed in sanitary container or as hos

INSTRUCTIONS FOR USE

REF

Product Ref. No.

LOT

Batch Number



Date of Manufacture



Use-by-date



Do Not Reuse



Sterilized using Ethylene Oxide and Single Sterile Barrier System



Caution



Consult Instructions for Use



Non-Pyrogenic



Do Not Use if package is Damaged and consult instructions for use



Do not resterilize



Keep away from sunlight



DEHP Free



Latex Free



Temperature limit

MD

Medical Device

UDI

Unique Device Identifier



Keep Dry



Non-Toxic/Non-Hazardous Substances



Importer



Distributer



Manufactured By:
Harsoria Healthcare Pvt. Ltd.
110 -111, Phase-IV, Udyog Vihar,
Gurgaon - 122 015, Haryana, INDIA.
Email: info@harsoria.com
Website: www.harsoria.com
Customer Care No.: +91 124 4523400


CE 0123

EU REP

Mdi Europa GmbH
Langenhagener Str. 71
30855 Langenhagen - GERMANY
Email: info@mdi-europa.com
Rev. Date :12.06.2025
Ref. Code: AW/IC_01a, Rev. 17


Color of Printing = Blue Pantone 2728cvc & Black
Color scheme = As per design
Base = Brown
Size = 698 mm X 320 X 178 mm
3 PLY, 1 kraft- 180 gsm (top Layer),
2 semi-kraft – 140 gsm

FILE NAME/TYPE	HEALFLON (TF_01) / FLON TYPE		 HARSORIA HEALTHCARE PVT. LTD. CMYK 	
DESIGNER	BALBEER SINGH		PANTONE COLORS USED BLACK BROWN	
MINIMUM	Current Size	Standard Size		
FONT TYPE AND SIZE	5.5	5		
	Arial	Time New Roman/Arial		
DATE STARTED	12/09/2025			
DATE COMPLETE	01/10/2025			
FILE LOCATION	Z:\CANNULA\FLON TYPE\HEALFLON (TF-01)\MC_01, Rev. 04.cdr		SIZE (LXWXH) - 698X320X178mm	



Color of Printing = Blue Pantone 2728cvc & Black
Color scheme = As per design
Base = Brown
Size = 625 mm X 275 X 204 mm
3 PLY, 1 kraft- 180 gsm (top Layer),
2 semi-kraft – 140 gsm
Only for 22-26 G

FILE NAME/TYPE	HEALFLON (TF_01) / FLON TYPE		 HARSORIA HEALTHCARE PVT. LTD. 			
DESIGNER	BALBEER SINGH		PANTONE			
MINIMUM	Current Size	Standard Size				
FONT TYPE AND SIZE	5.5	5	COLORS			
	Arial	Time New Roman/Arial				
DATE STARTED	12/09/2025		USED			
DATE COMPLETE	01/10/2025					
FILE LOCATION	Z:\CANNULA\FLON TYPE\HEALFLON (TF-01)\MC_01a, Rev. 03.cdr		SIZE (LXWXH) - 625X275X204mm			




Master Carton Sticker


Black print, Size= 100X125 mm

HEALFLON[®]

I.V. CANNULA WITH CATHETER
& INJECTION VALVE



USE IS RESTRICTED TO A QUALIFIED DOCTOR OR A TRAINED PARAMEDIC



STRIP STICKER



Manufactured By:

HARSORIA HEALTHCARE PVT. LTD.

110 -111, Phase-IV, Udyog Vihar,

Gurgaon - 122 015, Haryana, INDIA.

Email: info@harsoria.com

Website: www.harsoria.com

Customer Care No.: +91 124 4523400

M.L.No.: MFG/MD/2019/000135

EU REP

Mdi Europa GmbH

Langenhagener Str. 71

30855 Langenhagen - GERMANY

Email: info@mdi-europa.com

AW/MCS_01, Rev. 16

CE

0123

REF

LOT

UNITS

1000

Manufacturing Date

Expiry Date

UDI



(01)XXXXXXXXXXXXXXXXXX

(17)YYMMDD

(10)XXXXXX/XXXX

Variable Data Change as per Sales Order

FILE NAME/TYPE	HEALFLON (TF_01) / FLON TYPE		HARSORIA HEALTHCARE PVT. LTD. CMYK	
DESIGNER	BALBEER SINGH		PANTONE BLACK COLORS USED	
MINIMUM FONT TYPE AND SIZE	Current Size	Standard Size		
	5.5	5		
	Arial	Time New Roman/Arial		
DATE STARTED	12/09/2025			
DATE COMPLETE	01/10/2025			
FILE LOCATION	Z:\CANNULA\FLON TYPE\HEALFLON (TF-01)\MCS_01, Rev. 16.cdr		SIZE (LXH) - 100X125mm	

Actual Product Image

